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Engineering research approaches for hydrogen peroxide decontamination processes

藪田 啓奨* 杉山 弘和* 川崎 康司**
Keisho Yabuta Hirokazu Sugiyama Koji Kawasaki

This article presents the collaboration research between Airex and the University of Tokyo, which applies engineering research approaches for decontamination processes using hydrogen peroxide. Decontamination is a process to establish adequate environment in isolators or clean rooms in order to manufacture so-called sterile drug products such as injectables. Using hydrogen peroxide as a typical decontaminant, the process is critical regarding quality assurance of the produced products, but also productivity in manufacturing because the process requires long time duration. So far, experiments were performed using a simple isolator at Airex, and a statistical model was constructed based on the experimental results. The developed model can be used to propose the conditions of hydrogen peroxide loading and initial humidity that minimizes the total time of the decontamination process, i.e., sterilization and aeration. In future, the scope will be extended to cover more detailed aspects such as diffusion, reaction and condensation of hydrogen peroxide for improving the applicability of the model. The presented research would bring the following benefits (i) acceleration of isolator development for isolator makers, (ii) enhancement of productivity for pharmaceutical companies, and (iii) exploration of new research field for academia.

KEY WORDS: Pharmaceuticals, Isolator, Decontamination, Process modeling, Process design

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Fig. 1 Isolator for sterile manufacturing

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
Fig. 2 H₂O₂ concentration profile in isolator

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
Fig. 3 Relationship between design variables and objective functions considered in the model

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
Fig. 4 Contour plot of total time y_1 as an example of the model output¹⁰⁾

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	Current situation	This approach	Expected result
Isolator makers (e.g. Airex)	Experienced-based approach	Model development $y_i = f_i(x)$	Acceleration of isolator development
Pharmaceutical companies	Long time for decontamination		Enhancement of productivity
Academia	Needs for new research field		Exploration of new research field

Fig. 5 The role of this research in industry and academia

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